

Results

Figure 2 shows the TIC chromatogram. Samples A and B were made of the same polypropylene, but there were significant differences in the detected peaks.




Figure 2 TIC chromatograms

Figure 3 shows the difference analysis results of msFineAnalysis iQ. From the 11 peaks with an intensity ratio of 5% or more to the maximum peak, 2 peaks that were strong in sample A and 8 peaks that were strong in sample B were extracted.





Figure 3 Difference analysis result of msFineAnalysis iQ

Table 2 shows the peak list of the qualitative analysis results. The background color indicates the intensity difference between samples. Blue is strong in sample A, red is strong in sample B, and white is not different between samples. In Addition, to improve the accuracy of qualitative analysis, msFineAnalysis iQ searches for molecular ions from measured mass spectra. In this result, it was possible to confirm molecular ions in all 11 peaks.

Table 2 Peak list of qualitative analysis result

ID	General		Variance Component Analysis Result		Library Name	Total Result							Spectrum Info
	RT [min]	Height [%]	Class	Log2(B/A)		CAS#	Lib.	Similarity	Formula	DBE	MW	Molecular Weight Check	
001	1.33	63.19	A < B	1.54	n-Hexane	110-54-3	mainlib	877	C6 H14	0.0	86	✓	86
002	14.54	7.00	A = B	0.89	Pentadecane	629-62-9	mainlib	869	C15 H32	0.0	212	✓	212
003	14.74	18.25	A Only	< -4	2,4-Di-tert-butylphenol	96-76-4	mainlib	911	C14 H22 O	4.0	206	✓	206
004	14.77	69.39	A Only	< -4	Butylated Hydroxytoluene	128-37-0	mainlib	977	C15 H24 O	4.0	220	✓	220
005	15.65	9.66	B Only	> 4	Cetene	629-73-2	mainlib	857	C16 H32	1.0	224	✓	224
006	15.73	26.66	B Only	> 4	Hexadecane	544-76-3	mainlib	876	C16 H34	0.0	226	✓	226
007	16.85	27.70	B Only	> 4	Heptadecane	629-78-7	mainlib	865	C17 H36	0.0	240	✓	240
008	17.91	30.08	B Only	> 4	Octadecane	593-45-3	mainlib	822	C18 H38	0.0	254	✓	254
009	18.93	14.33	B Only	> 4	Nonadecane	629-92-5	mainlib	905	C19 H40	0.0	268	✓	268
010	19.19	100.00	A < B	2.94	7,9-Di-tert-butyl-1-oxaspiro(4,5)deca-6,9-diene-2,8-dione	82304-66-3	mainlib	924	C17 H24 O3	6.0	276	✓	276
011	19.89	7.99	B Only	> 4	Eicosane	112-95-8	mainlib	903	C20 H42	0.0	282	✓	282

Checked when molecular ion is detected

Figure 4 shows the mass spectra of ID004 (detected at 14.77 min in sample A) and ID010 (detected at 19.19 min in sample B). As a result of qualitative analysis, the former was a phenolic antioxidant Butylated hydroxytoluene (BHT). The latter was [7,9-Di-tert-butyl-1-oxaspiro(4,5)deca-6,9-diene-2,8-dione], which is a decomposition product of a hindered phenolic antioxidant Pentaerythritol tetrakis[3-(3',5'-di-tert-butyl-4'-hydroxyphenyl)propionate].




Figure 4 Mass spectra of ID004 and ID010

Conclusion

Using TD-100X and JMS-Q1600GC, it is possible to analyze additives in propylene products with high sensitivity. In addition, msFineAnalysis iQ can easily extract characteristic components from complex chromatograms and perform qualitative analysis with high accuracy. These devices and software are also expected to be useful in materials analysis.

Copyright © 2023 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

